

Pulsed NMR

Feb 28, 2017

Author: Jacob Cluff

Partners: Jason Pickering, Philip Rybak

Abstract:

Introduction:

Theory:

Procedure:

Results:

Figure 1:

Discussion:

Appendix:

1 Setup

2 Free-Induction Decay (FID)

A 90° Pulse FID

- i default oscilloscope settings: Time/Div 20 μs , CH1 1 V/cm, CH2 2V/cm, Trig: Ext, 1 V, POS.
- ii tune the oscilloscope until the FID (detector out, channel 1) is maximized.
- iii minimize the beat oscillations (mixer out, channel 2) by using the "Frequency Adjust".
- iv adjust "A-width" for 90° pulse, for maximum FID.
- v ensure optimum values by repeating steps (i), (ii), and (iii - this step).
- vi record all settings
- vii adjust the "Frequency Adjust" to see the beat oscillations in channel 2, sketch them, describe what causes them, and estimate T_2 - the time it takes for oscillation to reach 1/e of its original amplitude.

B Other Pulse Lengths

- i start with the "A-width" and "Frequency Adjust" settings for 90° pulse and beat oscillations from part A.
- ii increase "A-width" and describe what happens to the oscillation in channel 2.
- iii how is a 360° pulse width identified?
- iv sketch oscillations from channel 2 with pulse lengths of 90°, 180°, and 360° and explain the behavior.

3 Spin-Lattice Relaxation Time (T_1)

A Single Pulse

- i start with the 90° settings found in 2.A.v. Maximize "Detector Out" channel 1 by increasing "Repetition Time".
- ii record the "Repetition Time" that gives max amplitude for 90° pulse.
- iii decrease "Repetition Time" and explain reduction in amplitude.
- iv get an estimate for T_1 by the "Repetition Time" that gives about $1/e \approx 1/3$ of the maximum value.

B Two Pulses, Zero Crossing

- i Use the settings from 2.A.v for the in-tune 90° and set the repetition rate to about 100 ms.
- ii Turn the B pulses on (B pulses = 1) and start with a small delay - about 0.2 ms.
- iii Adjust the settings on the oscilloscope until both A and B pulses can be observed.
- iv Adjust the "A-width" to a 180° pulse (minimized)
- v Adjust the "B-width" to a 90° pulse (maximized)
- vi focus on the amplitude of the second pulse by switching "Sync" to B.
- vii find T_1 by adjusting the "Delay Time" (τ). T_1 is the time that minimized the amplitude.
- viii find a more statistically reliable value for T_1 ; measure the amplitude of the second pulse as a function of τ and fit the curve to a quadratic function and find the minimum of the curve.

4 Spin-Spin Relaxation Time (T_2)

A Two-Pulse Spin Echo

- i start with A and B pulses set as in 3.B

- ii adjust the "A-width" to maximized 90° pulse; the initial value of second amplitude should be minimized.
- iii adjust the "B-width" until it is maximized (180° pulse).
- iv get the second amplitude as a function of total delay time (2τ), using the cursor for Channel 1 to get the best values for the amplitude.
- v fit data to an exponential decay to find T_2 .

B Carr-Purcell Multiple Pulse Spin Echo Sequences

- i Start with same settings as 4.A, but start with 3 B pulses and make sure the CPMG switch is turned off.
- ii lower the repetition rate to avoid overheating of the rf amplifier.
- iii adjust the "Delay Time" and other possible oscilloscope settings until several time separated pulses can be observed.
- iv Add more B pulses and tweak the settings until a decay in their amplitudes can be seen.
- v obtain a function of the second amplitude as a function of "Delay Time" from the first 90° pulse).
- vi fit the data to an exponential decay to obtain a fitted value for T_2 .

C Meiboom-Gill Sequence

- i turn on the CPMG (Meiboom-Gill), repeat 4.B.iv - 4.B.vi, and obtain another value of T_2 ; why would the two numbers differ?

5 Find T_1 and T_2 using one of the experiments in a different medium (B glycerin or D petroleum jelly)